



# LRTM1 Antibody, Biotin conjugated

Product Code	CSB-PA884513LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBL6
Immunogen	Recombinant Human Leucine-rich repeat and transmembrane domain-containing protein 1 protein (1-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat and transmembrane domain-containing protein 1, LRTM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRTM1