




ZNF71 Antibody

Product Code	CSB-PA885697LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQZ8
Immunogen	Recombinant Human Endothelial zinc finger protein induced by tumor necrosis factor alpha protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endothelial zinc finger protein induced by tumor necrosis factor alpha (Zinc finger protein 71), ZNF71, EZFIT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF71

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA885697LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.