



LMBRD1 Antibody, HRP conjugated

Product Code	CSB-PA885730LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUN5
Immunogen	Recombinant Human Probable lysosomal cobalamin transporter protein (6-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable lysosomal cobalamin transporter (HDAg-L-interacting protein NESI) (LMBR1 domain-containing protein 1) (Nuclear export signal-interacting protein), LMBRD1, C6orf209 NESI
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	LMBRD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.