




ASNSD1 Antibody

Product Code	CSB-PA885745LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWL6
Immunogen	Recombinant Human Asparagine synthetase domain-containing protein 1 protein (258-357AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Asparagine synthetase domain-containing protein 1 (HCV NS3-transactivated protein 1), ASNSD1, NS3TP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ASNSD1

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA885745LA01HU at dilution of 1:100