



MRPL36 Antibody, Biotin conjugated

Product Code	CSB-PA885784LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0J6
Immunogen	Recombinant Human 39S ribosomal protein L36, mitochondrial protein (2-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L36, mitochondrial (L36mt) (MRP-L36) (BRCA1-interacting protein 1) (Mitochondrial large ribosomal subunit protein bL36m), MRPL36, BRIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPL36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.