



IFT80 Antibody, FITC conjugated

Product Code	CSB-PA885802LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2H3
Immunogen	Recombinant Human Intraflagellar transport protein 80 homolog protein (1-268AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 80 homolog (WD repeat-containing protein 56), IFT80, KIAA1374 WDR56
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IFT80
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.