



STOML1 Antibody, HRP conjugated

Product Code	CSB-PA887004LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBI4
Immunogen	Recombinant Human Stomatin-like protein 1 protein (79-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP), STOML1, SLP1 UNC24
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STOML1