



STOML1 Antibody, HRP conjugated

Product Code	CSB-PA887004LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBI4
Immunogen	Recombinant Human Stomatin-like protein 1 protein (79-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP), STOML1, SLP1 UNC24
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STOML1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.