



STOML1 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA887004LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UBI4 |
| Immunogen | Recombinant Human Stomatin-like protein 1 protein (79-290AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Stomatin-like protein 1 (SLP-1) (EPB72-like protein 1) (Protein unc-24 homolog) (Stomatin-related protein) (STORP), STOML1, SLP1 UNC24 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | STOML1 |