



TCEAL9 Antibody

Product Code	CSB-PA887044LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHQ7
Immunogen	Recombinant Human WW domain-binding protein 5 protein (1-104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription elongation factor A protein-like 9 (TCEA-like protein 9) (Transcription elongation factor S-II protein-like 9) (WW domain-binding protein 5) (WBP-5), TCEAL9, WBP5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TCEAL9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.