



RTP3 Antibody, Biotin conjugated

Product Code	CSB-PA887106LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQQ7
Immunogen	Recombinant Human Receptor-transporting protein 3 protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-transporting protein 3 (3CxxC-type zinc finger protein 3) (Transmembrane protein 7), RTP3, TMEM7 Z3CXXC3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTP3