



VPS25 Antibody, Biotin conjugated

Product Code	CSB-PA887108LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRG1
Immunogen	Recombinant Human Vacuolar protein-sorting-associated protein 25 protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein-sorting-associated protein 25 (hVps25) (Dermal papilla-derived protein 9) (ELL-associated protein of 20 kDa) (ESCRT-II complex subunit VPS25), VPS25, DERP9 EAP20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS25