



VPS25 Antibody, Biotin conjugated

Product Code	CSB-PA887108LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRG1
Immunogen	Recombinant Human Vacuolar protein-sorting-associated protein 25 protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein-sorting-associated protein 25 (hVps25) (Dermal papilla-derived protein 9) (ELL-associated protein of 20 kDa) (ESCRT-II complex subunit VPS25), VPS25, DERP9 EAP20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.