





TAPBPL Antibody

Product Code	CSB-PA887144LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BX59
Immunogen	Recombinant Human Tapasin-related protein (201-468AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tapasin-related protein (TAPASIN-R) (TAP-binding protein-like) (TAP-binding protein-related protein) (TAPBP-R) (Tapasin-like), TAPBPL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TAPBPL

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA887144LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA887144LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.