





SPZ1 Antibody

Product Code	CSB-PA887147LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXG8
Immunogen	Recombinant Human Spermatogenic leucine zipper protein 1 protein (181-430AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenic leucine zipper protein 1 (Testis-specific protein 1) (Testis-specific protein NYD-TSP1), SPZ1, TSP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SPZ1

Image



Western blot
 All lanes: SPZ1 antibody at 4µg/ml
 Lane 1: Mouse gonadal tissue
 Lane 2: NIH/3T3 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 50 kDa
 Observed band size: 50 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA887147LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.