



SPZ1 Antibody, HRP conjugated

Product Code	CSB-PA887147LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXG8
Immunogen	Recombinant Human Spermatogenic leucine zipper protein 1 protein (181-430AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenic leucine zipper protein 1 (Testis-specific protein 1) (Testis-specific protein NYD-TSP1), SPZ1, TSP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SPZ1