





RSPO3 Antibody

Product Code	CSB-PA887154LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXY4
Immunogen	Recombinant Human R-spondin-3 protein (22-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	R-spondin-3 (Protein with TSP type-1 repeat) (hPWTSR) (Roof plate-specific spondin-3) (hRspo3) (Thrombospondin type-1 domain-containing protein 2), RSPO3, PWTSR THSD2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	RSPO3

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA887154LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA887154LA01HU at dilution of 1:100