



RSPO3 Antibody, HRP conjugated

Product Code	CSB-PA887154LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXY4
Immunogen	Recombinant Human R-spondin-3 protein (22-272AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	R-spondin-3 (Protein with TSP type-1 repeat) (hPWTSR) (Roof plate-specific spondin-3) (hRspo3) (Thrombospondin type-1 domain-containing protein 2), RSPO3, PWTSR THSD2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	RSPO3