



RSPO3 Antibody, FITC conjugated

Product Code	CSB-PA887154LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXY4
Immunogen	Recombinant Human R-spondin-3 protein (22-272AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	R-spondin-3 (Protein with TSP type-1 repeat) (hPWTSR) (Roof plate-specific spondin-3) (hRspo3) (Thrombospondin type-1 domain-containing protein 2), RSPO3, PWTSR THSD2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	RSPO3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.