



NACA2 Antibody, Biotin conjugated

Product Code	CSB-PA887940LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H009
Immunogen	Recombinant Human Nascent polypeptide-associated complex subunit alpha-2 protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nascent polypeptide-associated complex subunit alpha-2 (Alpha-NAC-like) (Homo s 2.01) (Nascent polypeptide-associated complex subunit alpha-like) (NAC-alpha-like), NACA2, NACAL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NACA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.