



RPP21 Antibody, Biotin conjugated

Product Code	CSB-PA887986ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H633
Immunogen	Recombinant Human Ribonuclease P protein subunit p21 protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p21 (RNaseP protein p21) (EC 3.1.26.5) (Ribonuclease P/MRP 21 kDa subunit) (Ribonucleoprotein V), RPP21, C6orf135 CAT60
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPP21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.