



RPP21 Antibody, Biotin conjugated

Product Code	CSB-PA887986ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H633
Immunogen	Recombinant Human Ribonuclease P protein subunit p21 protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease P protein subunit p21 (RNaseP protein p21) (EC 3.1.26.5) (Ribonuclease P/MRP 21 kDa subunit) (Ribonucleoprotein V), RPP21, C6orf135 CAT60
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPP21