



EPS8L2 Antibody, HRP conjugated

Product Code	CSB-PA887991LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6S3
Immunogen	Recombinant Human Epidermal growth factor receptor kinase substrate 8-like protein 2 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epidermal growth factor receptor kinase substrate 8-like protein 2 (EPS8-like protein 2) (Epidermal growth factor receptor pathway substrate 8-related protein 2) (EPS8-related protein 2), EPS8L2, EPS8R2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EPS8L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.