





IL1RL2 Antibody

Product Code	CSB-PA888016LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB29
Immunogen	Recombinant Human Interleukin-1 receptor-like 2 protein (20-218AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 receptor-like 2 (IL-36 receptor) (IL-36R) (Interleukin-1 receptor-related protein 2) (IL-1Rrp2) (IL1R-rp2), IL1RL2, IL1RRP2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL1RL2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA888016LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA888016LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.