




EXOSC3 Antibody

Product Code	CSB-PA889083LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQT5
Immunogen	Recombinant Human Exosome complex component RRP40 protein (2-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exosome complex component RRP40 (Exosome component 3) (Ribosomal RNA-processing protein 40) (p10), EXOSC3, RRP40
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EXOSC3

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA889083LA01HU at dilution of 1:100