




EXOSC3 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA889083LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NQT5 |
| Immunogen | Recombinant Human Exosome complex component RRP40 protein (2-275AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Exosome complex component RRP40 (Exosome component 3) (Ribosomal RNA-processing protein 40) (p10), EXOSC3, RRP40 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | EXOSC3 |

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA889083LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.