



RBM22 Antibody, FITC conjugated

Product Code	CSB-PA889126LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NW64
Immunogen	Recombinant Human Pre-mRNA-splicing factor RBM22 protein (2-250AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-splicing factor RBM22 (RNA-binding motif protein 22) (Zinc finger CCCH domain-containing protein 16), RBM22, ZC3H16
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RBM22