





ELOVL2 Antibody

Product Code	CSB-PA889139LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXB9
Immunogen	Recombinant Human Elongation of very long chain fatty acids protein 2 protein (91-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 2 (EC 2.3.1.199) (3-keto acyl-CoA synthase ELOVL2) (ELOVL fatty acid elongase 2) (ELOVL FA elongase 2) (Very long chain 3-ketoacyl-CoA synthase 2) (Very long chain 3-oxoacyl-CoA synthase 2), ELOVL2, SSC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ELOVL2

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA889139LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA889139LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.