



PAQR5 Antibody, Biotin conjugated

Product Code	CSB-PA889141LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXK6
Immunogen	Recombinant Human Membrane progesterin receptor gamma protein (1-51AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane progesterin receptor gamma (mPR gamma) (Membrane progesterone P4 receptor gamma) (Membrane progesterone receptor gamma) (Progesterone and adipoQ receptor family member 5) (Progesterin and adipoQ receptor family member 5) (Progesterin and adipoQ receptor family member V), PAQR5, MPRG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PAQR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.