



TAS2R13 Antibody, HRP conjugated

Product Code	CSB-PA889153LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYV9
Immunogen	Recombinant Human Taste receptor type 2 member 13 protein (150-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 13 (T2R13) (Taste receptor family B member 3) (TRB3), TAS2R13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS2R13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.