



SPG21 Antibody, Biotin conjugated

Product Code	CSB-PA889158LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZD8
Immunogen	Recombinant Human Maspardin protein (1-308AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Maspardin (Acid cluster protein 33) (Spastic paraplegia 21 autosomal recessive Mast syndrome protein) (Spastic paraplegia 21 protein), SPG21, ACP33
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SPG21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.