



# HIGD1B Antibody, HRP conjugated

Product Code	CSB-PA889181LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P298
Immunogen	Recombinant Human HIG1 domain family member 1B protein (1-29AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HIG1 domain family member 1B (Protein CLST 11240), HIGD1B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HIGD1B