



TBC1D14 Antibody, Biotin conjugated

Product Code	CSB-PA889187LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2M4
Immunogen	Recombinant Human TBC1 domain family member 14 protein (136-435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 14, TBC1D14, KIAA1322
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.