



IMPACT Antibody, Biotin conjugated

Product Code	CSB-PA889191LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2X3
Immunogen	Recombinant Human Protein IMPACT protein (41-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein IMPACT (Imprinted and ancient gene protein homolog), IMPACT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	IMPACT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.