



PILRB Antibody, Biotin conjugated

Product Code	CSB-PA890735LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKJ0
Immunogen	Recombinant Human Paired immunoglobulin-like type 2 receptor beta protein (20-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired immunoglobulin-like type 2 receptor beta (Activating receptor PILR-beta) (Cell surface receptor FDFACT), PILRB, FDFACT
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PILRB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.