



# PRDM4 Antibody, HRP conjugated

Product Code	CSB-PA890920LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKN5
Immunogen	Recombinant Human PR domain zinc finger protein 4 protein (301-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PR domain zinc finger protein 4 (EC 2.1.1.-) (PR domain-containing protein 4), PRDM4, PFM1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRDM4