



PRDM4 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA890920LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UKN5 |
| Immunogen | Recombinant Human PR domain zinc finger protein 4 protein (301-500AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | PR domain zinc finger protein 4 (EC 2.1.1.-) (PR domain-containing protein 4), PRDM4, PFM1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PRDM4 |