



ZNRF3 Antibody, FITC conjugated

Product Code	CSB-PA890933HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULT6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase ZNRF3 protein (56-219AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase ZNRF3 (EC 2.3.2.27) (RING finger protein 203) (RING-type E3 ubiquitin transferase ZNRF3) (Zinc/RING finger protein 3), ZNRF3, KIAA1133 RNF203
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNRF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.