



TRAK1 Antibody, HRP conjugated

Product Code	CSB-PA891465LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPV9
Immunogen	Recombinant Human Trafficking kinesin-binding protein 1 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking kinesin-binding protein 1 (106 kDa O-GlcNAc transferase-interacting protein), TRAK1, KIAA1042 OIP106
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRAK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.