



PADI1 Antibody, HRP conjugated

Product Code	CSB-PA891543LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULC6
Immunogen	Recombinant Human Protein-arginine deiminase type-1 protein (1-663AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-arginine deiminase type-1 (EC 3.5.3.15) (Peptidylarginine deiminase I) (Protein-arginine deiminase type I), PADI1, PAD1 PDI1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PADI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.