



PADI3 Antibody, HRP conjugated

Product Code	CSB-PA891552LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW8
Immunogen	Recombinant Human Protein-arginine deiminase type-3 protein (1-664AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-arginine deiminase type-3 (EC 3.5.3.15) (Peptidylarginine deiminase III) (Protein-arginine deiminase type III), PADI3, PAD3 PDI3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PADI3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.