



# GABRQ Antibody, FITC conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA891564LC01HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                          |
| <b>Uniprot No.</b>         | Q9UN88   |
| <b>Immunogen</b>           | Recombinant Human Gamma-aminobutyric acid receptor subunit theta protein (348-611AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | FITC   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4       |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Gamma-aminobutyric acid receptor subunit theta (GABA(A) receptor subunit theta), GABRQ |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Others   |
| <b>Target Names</b>        | GABRQ  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.            |