




FBXW2 Antibody

Product Code	CSB-PA891730LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKT8
Immunogen	Recombinant Human F-box/WD repeat-containing protein 2 protein (72-340AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 2 (F-box and WD-40 domain-containing protein 2) (Protein MD6), FBXW2, FBW2 FWD2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW2

Image




Western Blot


Positive WB detected in: HepG2 whole cell lysate, K562 whole cell lysate, HeLa whole cell lysate, Mouse liver tissue, Mouse kidney tissue
All lanes: FBXW2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 52, 45 kDa
Observed band size: 52 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA891730LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA891730LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.