



FBXW2 Antibody, Biotin conjugated

Product Code	CSB-PA891730LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKT8
Immunogen	Recombinant Human F-box/WD repeat-containing protein 2 protein (72-340AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 2 (F-box and WD-40 domain-containing protein 2) (Protein MD6), FBXW2, FBW2 FWD2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.