



ACSL5 Antibody, Biotin conjugated

Product Code	CSB-PA891734HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULC5
Immunogen	Recombinant Human Long-chain-fatty-acid--CoA ligase 5 protein (33-683AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain-fatty-acid--CoA ligase 5 (EC 6.2.1.3) (Long-chain acyl-CoA synthetase 5) (LACS 5), ACSL5, ACS5, FACL5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSL5