



ACSL5 Antibody, HRP conjugated

Product Code	CSB-PA891734LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULC5
Immunogen	Recombinant Human Long-chain-fatty-acid--CoA ligase 5 protein (33-683AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain-fatty-acid--CoA ligase 5 (EC 6.2.1.3) (Long-chain acyl-CoA synthetase 5) (LACS 5), ACSL5, ACS5, FACL5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.