




HPCAL4 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA891785LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UM19 |
| Immunogen | Recombinant Human Hippocalcin-like protein 4 protein (2-191AA) |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:1000-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Hippocalcin-like protein 4 (HLP4), HPCAL4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | HPCAL4 |

Image



Western blot

All lanes: HPCAL4 antibody at 2µg/ml

Lane 1: Mouse brain tissue

Lane 2: Mouse kidney tissue

Lane 3: Mouse ovarian tissue

Lane 4: Human placenta tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 22 kDa

Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.