



pos5 Antibody

Product Code	CSB-PA891853HA01SXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UT98
Immunogen	Recombinant Schizosaccharomyces pombe NADH kinase pos5, mitochondrial protein (43-386AA)
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces pombe
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH kinase pos5, mitochondrial (EC 2.7.1.86), pos5
Immunogen Species	Schizosaccharomyces pombe
Research Area	Others
Target Names	pos5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.