



pos5 Antibody, Biotin conjugated

Product Code	CSB-PA891853HD01SXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UT98
Immunogen	Recombinant Schizosaccharomyces pombe NADH kinase pos5, mitochondrial protein (43-386AA)
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces pombe
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH kinase pos5, mitochondrial (EC 2.7.1.86), pos5
Immunogen Species	Schizosaccharomyces pombe
Research Area	Others
Target Names	pos5