



MAGED2 Antibody, Biotin conjugated

Product Code	CSB-PA891976LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNF1
Immunogen	Recombinant Human Melanoma-associated antigen D2 protein (38-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-associated antigen D2 (11B6) (Breast cancer-associated gene 1 protein) (BCG-1) (Hepatocellular carcinoma-associated protein JCL-1) (MAGE-D2 antigen), MAGED2, BCG1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAGED2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.